

Product datasheet for **TP500538**

Snrpd2 (NM_026943) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse small nuclear ribonucleoprotein D2 (Snrpd2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR200538 representing NM_026943 Red =Cloning site Green =Tags(s) MSLLNKPKSEMTPEELQKREEEEFNTGPLSVLTQSVKNNTQVLINCRNNKLLGRVKAFDRHCNMVLENV KEMWTEVPKSGKGGKKKSKPVNKDRYISKMFLRGDSVIVVLRNPLIAGK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	13.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_081219</u>
Locus ID:	107686
UniProt ID:	<u>P62317</u> , <u>Q14AF6</u>
RefSeq Size:	508
Cytogenetics:	7 A3
RefSeq ORF:	354



[View online »](#)

Synonyms: 1810009A06Rik; SMD2

Summary: Plays role in pre-mRNA splicing as core component of the SMN-Sm complex that mediates spliceosomal snRNP assembly and as component of the spliceosomal U1, U2, U4 and U5 small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome. Component of both the pre-catalytic spliceosome B complex and activated spliceosome C complexes. Is also a component of the minor U12 spliceosome.[UniProtKB/Swiss-Prot Function]