

Product datasheet for **TP500596**

6330403K07Rik (NM_134022) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RIKEN cDNA 6330403K07 gene (6330403K07Rik), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR200596 protein sequence Red =Cloning site Green =Tags(s)
	MAKGPFSRKPWGPEQMYFDS DPEYEG LFDKPPLEEGHTARAPKSASSAGRKSGRRSSGRAPGTRAGLSRK AAVRPEPK EEEPPVEEGCYLDHFP HLSIFIYAAIAFSITSCIFTYIHLQLA
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	13.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_598783</u>
Locus ID:	103712
UniProt ID:	<u>Q8VBZ2</u> , <u>Q6NSQ3</u>
RefSeq Size:	1626
Cytogenetics:	11 B4
RefSeq ORF:	366



[View online »](#)

Synonyms: AI415455; AV161034; UGS148