

## Product datasheet for **TP500608**

### InsI3 (NM\_013564) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse insulin-like 3 (InsI3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR200608 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MRAPLLLMLLALGSALRSPQPPEARAKLCGHHLVRTLVRVCGGPRWSPEATQPVETRDRELLQWLEQRHL LHALVADVDPALDPQLPLQASQRQRSAATNAVHRCCLTGCTQQDLLGLCPH
	<b>TR</b> TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	13.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_038592</a>
Locus ID:	16336
UniProt ID:	<a href="#">Q5RL10</a>
RefSeq Size:	587
Cytogenetics:	8 34.43 cM
RefSeq ORF:	369



[View online »](#)

Synonyms: Rlf; Rlnl