

## Product datasheet for **TP500620**

### Vpreb3 (NM\_009514) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse pre-B lymphocyte gene 3 (Vpreb3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from Mouse cDNA ORF Clone, MR200620, encoding Mouse full-length Vpreb3.
Tag:	C-MYC/DDK
Predicted MW:	13.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_033540</a>
Locus ID:	22364
UniProt ID:	<a href="#">Q61243</a>
RefSeq Size:	602
Cytogenetics:	10 C1
RefSeq ORF:	369
Synonyms:	8HS-20; AI528709; Vpreb-3



[View online »](#)