

Product datasheet for TP500842

Ndufb7 (NM_025843) Mouse Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Mouse NADH:ubiquinone oxidoreductase subunit B7 (Ndufb7), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug

Species: Mouse

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >MR200842 protein sequence
Red=Cloning site Green=Tags(s)

MG AHLTRRYLWDASVEPDPEKIPSFPPDLGFPERKERV MVATQQEMMDAQLTLQQRDYCAHYLIRLLKCK
RDSFPNFLACKHEQHDWDYCEHLDYVKRMKEFERERLLQRKKRRALKEARVAQGQGEVGP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-MYC/DDK

Predicted MW: 16.3 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C after receiving vials.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_080119](#)

Locus ID: 66916

UniProt ID: [Q9CR61](#)

RefSeq Size: 533

Cytogenetics: 8 C2

RefSeq ORF: 414



[View online »](#)

Synonyms: 1110002H15Rik

Summary: Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.[UniProtKB/Swiss-Prot Function]