

## **Product datasheet for TP500888**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## 5730437N04Rik (BC019420) Mouse Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Mouse RIKEN cDNA 5730437N04 gene (cDNA clone

MGC:30293 IMAGE:5053733), complete cds, with C-terminal MYC/DDK tag, expressed in

HEK293T cells, 20ug

**Species:** Mouse

**Expression Host:** HEK293T

Expression cDNA Clone or AA Sequence:

>MR200888 protein sequence Red=Cloning site Green=Tags(s)

METSGPGPGESSELEAPGSPDDRLFLVKGGIFLGSAAAAGMLAGFVTTLSLAKKKSPEWFNKGTMATAAL PESGSSLALRALGWGSLYAWCGVGVISFAVWKALGVHSMKDFRSKMQSIFPPIPKNHESAEEWEEVLKWK

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV** 

Tag: C-MYC/DDK

Predicted MW: 14.9 kDa

**Concentration:** >0.05 μg/μL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

**Storage:** Store at -80°C after receiving vials.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

Locus ID: 70544
UniProt ID: Q8VCR3

RefSeq Size: 939 Cytogenetics: 17 A1

RefSeq ORF: 420







**Synonyms:** 1110008A10Rik, 2310046K16Rik