

Product datasheet for **TP501071**

Gtsf1l (NM_026630) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse gametocyte specific factor 1-like (Gtsf1l), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from Mouse cDNA ORF Clone, MR201071, encoding Mouse full-length Gtsf1l.
Tag:	C-MYC/DDK
Predicted MW:	17.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_080906
Locus ID:	68236
UniProt ID:	Q9CWD0
RefSeq Size:	880
Cytogenetics:	2 H2
RefSeq ORF:	453
Synonyms:	2410116G06Rik



[View online »](#)