

Product datasheet for **TP501146**

Hspb3 (NM_019960) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse heat shock protein 3 (Hspb3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201146 protein sequence Red =Cloning site Green =Tags(s)
	 MAKIILRHLIETPVRYQEEFEARGLEDLCRLDHTLYALPGPTIEDLSKARGAGTPQALAEDSASTEKPPGE GKSRFQILLDVVQFLPEDIIIQTFEGWLLIKAQHGTRMDEHGFISRSFTRQYKLPDGVETKDL SAILCHD GILVVEVKDSLGTK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	17.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_064344
Locus ID:	56534
UniProt ID:	Q9QZ57
RefSeq Size:	781
Cytogenetics:	13 D2.2



[View online »](#)

RefSeq ORF: 462

Synonyms: 2310035K17Rik; AI844863; Hsbp3; spb3

Summary: Inhibitor of actin polymerization.[UniProtKB/Swiss-Prot Function]