

## Product datasheet for **TP501207**

### **Btf3l4 (NM\_027453) Mouse Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Mouse basic transcription factor 3-like 4 (Btf3l4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence from Mouse cDNA ORF Clone, MR201207, encoding Mouse full-length Btf3l4.
<b>Tag:</b>	C-MYC/DDK
<b>Predicted MW:</b>	17.7 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C after receiving vials.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_081729</a>
<b>Locus ID:</b>	70533
<b>UniProt ID:</b>	<a href="#">Q9CQH7</a> , <a href="#">A2A7Z4</a>
<b>RefSeq Size:</b>	3111
<b>Cytogenetics:</b>	4 C7
<b>RefSeq ORF:</b>	474
<b>Synonyms:</b>	4632412E09Rik; 5730434I03Rik; AI316863; AU022207; AW552446



[View online »](#)