

Product datasheet for TP501214

Trappc6a (NM_025960) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse trafficking protein particle complex 6A (Trappc6a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201214 protein sequence Red =Cloning site Green =Tags(s) MADAVLFEFLHTEMVAELWAPDPDPGSGGKRRSLSVLEGLGFRVGQALGERLPLETPAFREELDALKFLC RDLWAAMFQKHMDGLRTNHQGTIVLQDNSFPLLVTMGSGPQYLEEAPKFLAFTCGLLCGALHTLGFQS LV TASVASLPACKFQVVIQKS TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	17.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_080236</u>
Locus ID:	67091
UniProt ID:	<u>Q78XR0</u>
RefSeq Size:	782


[View online »](#)

Cytogenetics:	7 A3
RefSeq ORF:	477
Synonyms:	1810073E21Rik; 4930519D19Rik; AI480686; mhyp; TRS33
Summary:	May play a role in vesicular transport during the biogenesis of melanosomes. [UniProtKB/Swiss-Prot Function]