

Product datasheet for **TP501338**

Snx12 (NM_018875) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse sorting nexin 12 (Snx12), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201338 protein sequence Red =Cloning site Green =Tags(s)
	 MSDTAVADTRRLNSKPQDLTDAYGPPSNFLEIDIFNPQTVGVGRARFTTYEVRMRTNLPFKLKESCVR RYSDFEWLKNELERDSKIVVPLPGEALKRQLPFRGDEGIFEESFIEERRQGLEQFINKIAGHPLAQNER CLHMFLQEEAIDRNYVPGKVLGEKDC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	19.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_061363
Locus ID:	55988
UniProt ID:	O70493 , Q3TGS7
RefSeq Size:	1000
Cytogenetics:	X D



[View online »](#)

RefSeq ORF: 501

Synonyms: 2610001F05Rik; AW045757; SDP8

Summary: May be involved in several stages of intracellular trafficking.[UniProtKB/Swiss-Prot Function]