

Product datasheet for TP501365

Cetn3 (NM_007684) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse centrin 3 (Cetn3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201365 protein sequence Red=Cloning site Green=Tags(s)
	MSLALRGELVVDKTKRKKRRELSEEQKQEIKDAFELFDTDKDQAIDYHELKVAMRALGFDVKKADVLKIL KDYDREATGKITFEDFNQVTDWILERDPHEEILKAFKLFDDDDSGKISLRNLRVARELGENMSDEELR AMIEEFDKGDGEINQEEFIAIMTGD
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	19.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_031710
Locus ID:	12626
UniProt ID:	O35648 , Q545L8
RefSeq Size:	1013
Cytogenetics:	13 C3



[View online »](#)

RefSeq ORF: 504

Synonyms: MmCEN3

Summary: Plays a fundamental role in microtubule-organizing center structure and function.
[UniProtKB/Swiss-Prot Function]