

Product datasheet for TP501442

Ttc9c (NM_027412) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse tetratricopeptide repeat domain 9C (Ttc9c), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201442 protein sequence Red =Cloning site Green =Tags(s) MEKRLQEAQVYKEEGNQRYREGKYRDAVSRYHRALLQLRGLDPSLPSPLSSLGPQGPALTPEQENILHTI QTHCYNNLAACLLQMEPVNYERVREYSQKVLERQPDNAKALYRAGVAFFHLQDYDRARHLLAAVNRP PK DANVRRYLQLTQSELSSYHRKEKQLYLGMFG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	20 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_081688
Locus ID:	70387
UniProt ID:	Q810A3
RefSeq Size:	3470


[View online »](#)

Cytogenetics: 19 A

RefSeq ORF: 513

Synonyms: 2210019E14Rik; 6330408J23Rik