

## Product datasheet for **TP501448**

### **C530044N13Rik (BC034347) Mouse Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RIKEN cDNA C530044N13 gene (cDNA clone MGC:27760 IMAGE:2654874), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201448 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MERSILLRPRRRHSVWTDEGLVLNSQFLYDASRCPALKQAQDHWLDQQLNIAEQKQCQHAIVFQHIPLFL QSIDEDDDYFNLTkTVRkELAEKLTRAGIRAVFSGHYHRNAGGTYQNLDMMVSSAIGCQLGKDTHTGLRVV AITAEKIVHRYYSLELSQGGVEEDLKELLKE
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	19.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	223978
UniProt ID:	<a href="#">Q8BFS6</a>
RefSeq Size:	2343
Cytogenetics:	16 A1



[View online »](#)

RefSeq ORF: 516

Synonyms: MGC27760

Summary: Protein phosphatase that dephosphorylates AKT family kinase specifically at 'Ser-473', blocking cell cycle progression and promoting cell apoptosis. May play an inhibitory role in glucose uptake by adipocytes (By similarity).[UniProtKB/Swiss-Prot Function]