

Product datasheet for TP501465

2310061I04Rik (BC100355) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RIKEN cDNA 2310061I04 gene (cDNA clone MGC:117734 IMAGE:30527188), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201465 protein sequence Red =Cloning site Green =Tags(s) MEEHLAVMHERLRQELPTLFLRSHDYTIYSMDVEFINEILNIRTKGRTFYVMSLTLCRFLVWNYFAQFRL EILQLTRHPENWTLQARWRLIGLPIHMLFLRFYRRDKEDLYRTFADFSTFYLNSSGLICRHHLDKLMPSH SPSTPVKLLLVGALVALGLSEPEPSLHLC SKA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	20.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	69662
RefSeq Size:	893
Cytogenetics:	17 B1
RefSeq ORF:	516



[View online »](#)