

Product datasheet for TP501728

Ift22 (NM_026073) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse intraflagellar transport 22 (Ift22), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201728 protein sequence Red =Cloning site Green =Tags(s)
	 MLKAKILFVGPCESGKTVLNFLTSSDITEYNPTQGVRIEFENPHVTSNNKGTGCEFELWDCGGDSKF ESCWPALMKDAHGWWVFNADIPSHLKEIEMWYSCFVQQFLQDSHCMLVAHHKPGSGGERGSLALSPPL NKLKLVHSNLEEDPEEVRVEFIKYLKSIINSMSERDREEMLIIT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	20.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_080349</u>
Locus ID:	67286
UniProt ID:	<u>Q9DAI2</u>
RefSeq Size:	911
Cytogenetics:	5 G2



[View online »](#)

RefSeq ORF: 558

Synonyms: 3110017O03Rik; AI835745; AI847427; Rabl5

Summary: Small GTPase-like component of the intraflagellar transport (IFT) complex B.[UniProtKB/Swiss-Prot Function]