

Product datasheet for TP501729

Ndufaf3 (NM_023247) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse NADH:ubiquinone oxidoreductase complex assembly factor 3 (Ndufaf3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201729 protein sequence Red =Cloning site Green =Tags(s)

MATALGFRCLFRTRPAPLCRHVDRLWRNPRRGHRLSPADDELYQRTRISLLQSEFPQAVVIDSYSSRGFT
ICGNRVFGPCVLLPQTVVHWNVGS HQDITEESFSLFWMLEPRIEIVVGTGNKTERLHPQVLQAMRQRGI
AVEVQDTPNACATFNFLCHEGRVTGAALIPPPGETALASLGQASE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	20.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_075736</u>
Locus ID:	66706
UniProt ID:	<u>Q9JKL4</u>
RefSeq Size:	1092
Cytogenetics:	9 F2



[View online »](#)

RefSeq ORF: 558

Synonyms: 4733401H18Rik

Summary: Essential factor for the assembly of mitochondrial NADH:ubiquinone oxidoreductase complex (complex I).[UniProtKB/Swiss-Prot Function]