

Product datasheet for TP501736

Tex29 (NM_029326) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse testis expressed 29 (Tex29), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201736 protein sequence Red=Cloning site Green=Tags(s)

MKDTKEIKRSPPHLLKKFAVCDIPLYDICDYNVTRERCRLDCCFYRGVCYEKAVPIYVQVFFTLIWFVA
GAFIIAVIYRVIQGTKKEKKPLLPEDSQVEILKSPTPELIPEPIPEPEPIREPPPPVKKTESPI
EAGCCLWMKSKPAKDQPQKETAPEPPSNPEVKKVNSGSAVPQAAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	20.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_083602
Locus ID:	75528
UniProt ID:	Q9DA77
RefSeq Size:	805
Cytogenetics:	8 A1.1



[View online »](#)

RefSeq ORF: 561

Synonyms: 1700018L24Rik