

Product datasheet for TP501771

1810063B07Rik (BC027644) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RIKEN cDNA 1810063B07 gene (cDNA clone MGC:41347 IMAGE:1248085), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201771 protein sequence Red =Cloning site Green =Tags(s)
	MEQRLAEFREARKRASLVAQPSTSSQSVQTSGAKAEPAAATPKTATGWLTRFLKRKANPAIAQAQPNQPQ EAGQQLPESTAVPLPSSCRQSFLTNI FLK VLLWLVLLGLFVELEFGLAYFVLSMFYWMYVGTRGPEEK EGEKSAYSVFNPGCEAIQGTLTAEQLEQELQLRPPQGSRTSPSCSSYP
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	20.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	67509
UniProt ID:	Q8K190
RefSeq Size:	758
Cytogenetics:	14 A3



[View online »](#)

RefSeq ORF: 564

Synonyms: 1810063B07Rik; 4930488P03Rik