

Product datasheet for TP501798

Stmn4 (NM_019675) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse stathmin-like 4 (Stmn4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201798 representing NM_019675 Red =Cloning site Green =Tags(s)
	 MTLAAYKEKMKELPLVSLFCSCFLSDPLNKSSYKYEADTVDLNWCVISDMEVIELNKCTSGQSFEVILKP PSFDGVPFNASLPRRRDPSLEEIQKKLEAAEERRKYQEAELLKHLAEKREHEREVIQKAIEENNNFIKM AKEKLAQKMESNKENREAHLAAMLRLQEKDKHAEVVRKNKELKEEASR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	22.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_062649
Locus ID:	56471
UniProt ID:	P63042
RefSeq Size:	1328
Cytogenetics:	14 D1



[View online »](#)

RefSeq ORF: 567

Synonyms: Lagl; R75359; RB3

Summary: Exhibits microtubule-destabilizing activity.[UniProtKB/Swiss-Prot Function]