

Product datasheet for **TP501865**

Myl4 (NM_010858) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse myosin, light polypeptide 4 (Myl4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201865 representing NM_010858 Red =Cloning site Green =Tags(s) MPPKKPEPKKETAKPAAAPAPAASAAPEPLKDSAFDPKSVKIDFSADQIEEFKEAFSLFDRTPTGEMKIT YGQCGDVLRLALGQNPTNAEVLRLVGLGPKPEEMSSKTLDFEMFLPILQHISRNKEQGTIEDFVEGLRVFDK ESNGTVMGAELRHVLTATLGEKMSEAEVEQLLSGQEDANGCINYEAFVKHIMSG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	21.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_034988</u>
Locus ID:	17896
UniProt ID:	<u>Q9CZ19</u>
RefSeq Size:	883
Cytogenetics:	11 67.79 cM



[View online »](#)

RefSeq ORF: 579

Synonyms: ALC1; AMLC; ELC; ELC1a; GT1; MLC1a; Myla