

Product datasheet for TP501995

Med29 (NM_026042) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse mediator complex subunit 29 (Med29), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR201995 protein sequence Red =Cloning site Green =Tags(s) MAAPQPQAAAVSSASGVSGPGSAGGPGPQQPQPTQLVGSAQSGLLQQQQQDFDPVQRYKMLIPQLK ESL QTLMKVAAQNLIQNTNIDNGQKSSDAPLQRFDKCLEEFYALCDQLELCLRLAHECLSQSCDSAKHSPTLV PTATKPDVQPDSPYPQYLAVIKAQITCAKDIHTALLDCANKVTGKTTAPSTGPGGSL TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	21 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_080318
Locus ID:	67224
UniProt ID:	Q9DB91
RefSeq Size:	2162


[View online »](#)

Cytogenetics:	7 A3
RefSeq ORF:	597
Synonyms:	2810405O22Rik; AU020902; Ixl
Summary:	<p>Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors (By similarity). [UniProtKB/Swiss-Prot Function]</p>