

Product datasheet for TP502006

Dnajc5b (NM_001163537) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse Dnaj heat shock protein family (Hsp40) member C5 beta (Dnajc5b), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR202006 representing NM_001163537 Red=Cloning site Green=Tags(s)
	MACNAPNQRQRTLSTSGESLYEILGLHKGASCEEIKKTYRKLALRHHPDKNPDDPSAAEFKEINNAHTI LTDTSKRNIYDKYGLGLYVAEQFGDENVNTYFMLSSWWAKTLFIIIGLLTGCYFCCCLCCCCNCCCGHC RPKSSTPEEEFYVSPEDLEEQIRTDMEKDMDFPVVLQPTNTNEKTQLIREGSRSYCTDS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	23.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001157009</u>
Locus ID:	66326
UniProt ID:	<u>Q9CQ94</u>
RefSeq Size:	960
Cytogenetics:	3 A2



[View online »](#)

RefSeq ORF: 597

Synonyms: 1700008A05Rik