

Product datasheet for **TP502013**

TagIn3 (NM_019754) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse transgelin 3 (TagIn3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR202013 protein sequence Red =Cloning site Green =Tags(s)
	<p>MANRGPSYGLSREVQEKIEQKYDADLENKLVDWIILQCAEDIEHPPPGRAHFQKWLMDGTVLCKLINSLY PPGQEPKISESKMAFKQMEQISQFLKAAEVYGVRTTDIFQTVDLWEGKDMAAVQRTLALGSAVATKD DGCYRGEPSWFHRKAQQNRRGFSEEQLRQGQNVIGLQMGSNKGASQAGMTGYGMPRQIM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	22.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_062728</u>
Locus ID:	56370
UniProt ID:	<u>Q9R1Q8</u>
RefSeq Size:	1269
Cytogenetics:	16 B5



[View online »](#)

RefSeq ORF: 600

Synonyms: 2700038H05Rik; 2900005O10Rik; A1426007; Np25