

## Product datasheet for TP502092

### Commd10 (NM\_178377) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse COMM domain containing 10 (Commd10), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR202092 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MAASAALILPESPSMKKAVPLINAIDTGRFPRLSRILQKLHLKAESSFSEEEEEKLQAAFSLEKQELHL VLETISFVLEQAVYHNVKPAALQQQLEMIHLRKDKAEAFASAWSAMGQETVEKFRQRILGPHKLETVGWQ LNLQMAHSAQAKLQSPQAVLQLGVSKEDAKNVEKVLVEFNHKELFDFYNKLETIQAQLDSL  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	22.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_848464</a>
Locus ID:	69456
UniProt ID:	<a href="#">Q8JZY2</a> , <a href="#">Q3TKN7</a>
RefSeq Size:	1550
Cytogenetics:	18 C



[View online »](#)

RefSeq ORF: 609

Synonyms: 2310003A05Rik; AU019438; DRWMS2

Summary: May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes. May down-regulate activation of NF-kappa-B.[UniProtKB/Swiss-Prot Function]