

## Product datasheet for **TP502267**

### Abhd14b (NM\_029631) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse abhydrolase domain containing 14b (Abhd14b), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR202267 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MAGVDQHEGTIQVQGQNLFFRETRPGSGQPVRFSVLLLHGIRFSSETWQNLGTLQRLAEAGYRAVAIDLPLGLGRSKEAAAPAPIGEPAPGSFLAAVDTLELGPPVVISPSLSGMYSLPFLVAPGSQLRGFVVPVAPICTDKINAVDYASVKTPALIVYGDQDPMGSSSFQHLKQLPNHRVLMGAGHPCYLDKPDEWHKGLLDLFLQGLA

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-MYC/DDK
Predicted MW:	22.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_083907</a></u>
Locus ID:	76491
UniProt ID:	<u><a href="#">E9QN99</a></u>
RefSeq Size:	1490
Cytogenetics:	9 F1



[View online »](#)

RefSeq ORF: 633

Synonyms: 1810013B01Rik