

Product datasheet for **TP502371**

Adamts15 (BC034843) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse ADAMTS-like 5 (cDNA clone MGC:41037 IMAGE:2938098), complete cds, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR202371 representing BC034843 Red =Cloning site Green =Tags(s)

MAIHPDPALHCLLVTLCSLKPRLFPLLPFLLWTLTLLSCLDLRGAAQGPGEWTPWGSWSRCSSSCGRGVSVRR
RQCVRLPGEELCWGDSHEYRVCQLPDCPPGTTTFRDLQCSLYNGHPVLDTQKTYQWVPFHGAPNVCDLNC
LAEGHAFYHSFGRVLDGTPCTPGTQGLCVAGRCLSAGCDGILGSGTLEDRCGLCGGANDSCLFVQRVFRD
AAIGP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	76.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
Locus ID:	66548
RefSeq Size:	2099
Cytogenetics:	10 C1
RefSeq ORF:	645



[View online »](#)

Synonyms: MGC41037