

## Product datasheet for **TP502504**

### Vstm2l (NM\_198627) Mouse Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Mouse V-set and transmembrane domain containing 2-like (Vstm2l), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug

**Species:** Mouse

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >MR202504 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MGAPLAAALGALHYLALFLQLGGATRPAGHAPWDNHVSGHALFTETPHDMMTARTGEDVEMACSFRRGSGSP  
SYSLEIQWWYLRSHRDWTDKQTWASNQLKASQQEDSGKDATKISVVKVGSNISHKLRLSRVKPTDEGTY  
ECRVIDFSDGGRGVPRVLCLLIPLPAPPAPRPRGQPPGEEPGRGPTLLFLIILPGTGSSTPREAEPHQ  
HAGGCPARQS

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-MYC/DDK

**Predicted MW:** 23.7 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C after receiving vials.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_941029](#)

**Locus ID:** 277432

**UniProt ID:** [Q6PDS0](#), [A2ACD2](#)

**RefSeq Size:** 1381



[View online »](#)

Cytogenetics: 2 H1

RefSeq ORF: 663

Synonyms: Gm691