

Product datasheet for TP502759

BC005537 (NM_024473) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse cDNA sequence BC005537 (BC005537), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR202759 protein sequence Red =Cloning site Green =Tags(s) MGDPNRKKQALNRLRAQLRKKKESLADQDFDKMYIAFVFKEKKKKKSALFEVSEVIPVMTNNYEENILKG VRDSSYSLESSLELLQKDVVQLHAPRYQSMRRDVIGCTQEMDFILWPRNDIEKIVCLLFSRWKESDEPFR PVQAKFEFHGHDYEKQFLHVLRSRKDKTGIVVNNPNQSVFLFIDRQHLQTPKNKATIFKLCSICLYLPQEQ LTHWAVGTIEDHLRPYMPE TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	27.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_077793
Locus ID:	79555
UniProt ID:	Q99LU8
RefSeq Size:	2551


[View online »](#)

Cytogenetics: 13 A3.1

RefSeq ORF: 687

Synonyms: 8030460C05Rik; C78201