

## Product datasheet for TP502879

### D630039A03Rik (NM\_178727) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RIKEN cDNA D630039A03 gene (D630039A03Rik), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR202879 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MRMKGLACPCPALPKFLELGTCLMAEGSRTQAPGQPPVNIQFLRAQYEGLRKQQRQAHLMVFPKEG TM PTPAESMTSAVWINKERRSSLSP EETDSEAEGMLEEGDRSCHQASETPWHTYLELHRLVQTFHLEAGHQG NPKGYLVEPEPRFSPEGDTDVLKNNQKTQQEVKIPEVAQCQSQEGCAPLQVAGSRAQYLSSTKHLRSGKP AYYFPQRKPPRISQAARNLGLYGPP
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	26.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_848842</a></u>
Locus ID:	242484
UniProt ID:	<u><a href="#">Q8K0M7</a></u>



[View online »](#)

RefSeq Size: 2929  
Cytogenetics: 4 B3  
RefSeq ORF: 708  
Synonyms: RP23-334A5.4