

## Product datasheet for **TP503027**

### **Snrnp35 (NM\_029532) Mouse Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Mouse small nuclear ribonucleoprotein 35 (U11/U12) (Snrnp35), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	A DNA sequence from Mouse cDNA ORF Clone, MR203027, encoding Mouse full-length Snrnp35.
<b>Tag:</b>	C-MYC/DDK
<b>Predicted MW:</b>	29.7 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C after receiving vials.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_083808</a>
<b>Locus ID:</b>	76167
<b>UniProt ID:</b>	<a href="#">Q9D384</a>
<b>RefSeq Size:</b>	1133
<b>Cytogenetics:</b>	5 F
<b>RefSeq ORF:</b>	732
<b>Synonyms:</b>	6330548G22Rik



[View online »](#)