

## Product datasheet for **TP503063**

### Tspan6 (NM\_019656) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse tetraspanin 6 (Tspan6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203063 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MASPSRRLQTKPVITCLKSVLLIYTFIFWITGVILLAVGIWGVKVSLENYFSLNEKATNVPFVLIGTGTV IILLGTFGCFATCRASAWMLKLYAMFLTLIFLVELVAAIVGFVFRHEIKNSFKSNYENALKEYNSTGDYR SEAVDKIQSTLHCCGVTNYGDWKGTYNYSETGFPKSCCKLEGYCPQRDADKVNEEGCFIKVMTTIESEM VWAGISFGVACFQLIGIFLAYCLSRAITNNQYEIV  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	27.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_062630</a>
Locus ID:	56496
UniProt ID:	<a href="#">Q99L96</a>
RefSeq Size:	1763



[View online »](#)

Cytogenetics: X E3

RefSeq ORF: 738

Synonyms: 6720473L21Rik; AI316786; T245; Tm4sf; Tm4sf6; TSPAN-6