

Product datasheet for **TP503151**

Ankrd45 (NM_028664) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse ankyrin repeat domain 45 (Ankrd45), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203151 protein sequence Red =Cloning site Green =Tags(s)

MEPEETLESESEKSLFSSQQEYEESQEAETGAENPLLQPTLTGDVEGLQKIFEDPEHPHHEHAVQLLL
EEDIVGRNLLYAACMAGKSDVIKALAKYGVNLNEATARGYTLLHCAAAGRLETALKALVELDVDIEALNF
RGEKARDVAARYSQVECVNFLDWADARLILKKIITKSSLIITDPEKGPGKLFKEDKSTILNACRLKNEWL
ESHPEASISEIFEQKQLEDIVSPILAKMSTPRHFAAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	27.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_082940</u>
Locus ID:	73844
UniProt ID:	<u>Q810N6</u>
RefSeq Size:	3723



[View online »](#)

Cytogenetics:	1 H2.1
RefSeq ORF:	744
Synonyms:	4933409K03Rik