

## Product datasheet for TP503212

### Hdhd3 (NM\_024257) Mouse Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Mouse haloacid dehalogenase-like hydrolase domain containing 3 (Hdhd3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug

**Species:** Mouse

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >MR203212 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MAHRLQMRLLTWDVKDTLIKLRPVGEEYASKARAHGVVEDITVEQAFRQAYRAQSHNFPNYGLSRGLT  
SRQWWKDVVLTFRLAGVDPDAQAMTPVADQLYEDFSSPFTWQVLEGAEMTLKGCRRGLKLAVVSNFD  
RR  
LEDILTGLGLREHFDFVLTSEAVGCPKPDPRIFREALQRACVEPAVAHVGDSYLCDYQGSQAVGMHSFL  
VAGSEPLDSAVRDSVPKEHILPSLSHLLPALDLLEASSPMS

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-MYC/DDK

**Predicted MW:** 28 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C after receiving vials.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_077219](#)

**Locus ID:** 72748

**UniProt ID:** [Q9CYW4](#)



[View online »](#)

RefSeq Size: 1098  
Cytogenetics: 4 B3  
RefSeq ORF: 753  
Synonyms: 2810435D12Rik