

Product datasheet for **TP503289**

H2-Aa (NM_010378) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse histocompatibility 2, class II antigen A, alpha (H2-Aa), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203289 protein sequence Red =Cloning site Green =Tags(s)
	<p>MPRSRALILGVLALTTMLSLCGGEDDIEADHVGTYGISVYQSPGDIGQYTFEFDGDELFFVVDLDDKKEVW MLPEFGQLASFDPPQGGLQNIQAVVVKHNLGVLTKRSNSTPATNEAPQATVFPKSPVLLGQPNTLICFVDNIF PPVINITWLRNSKSVADGVYETSEFFVNRDYSFHKLSYLTIFPSDDDIYDCKVEHWGLEEPVLKHWEPEIP APMSELTETVVCALGLSVGLVGIIVGTIFIIQGLRSGGTSRHPGPL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	28.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_034508
Locus ID:	14960
UniProt ID:	P14434
RefSeq Size:	1084



[View online »](#)

Cytogenetics: 17 17.98 cM

RefSeq ORF: 771

Synonyms: Aalpha; H-2Aa; H2Aa; I-Aalpha; Ia-1; Ia1; IAalpha