

## Product datasheet for TP503336

### Fam221a (NM\_001172216) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse family with sequence similarity 221, member A (Fam221a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203336 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MERLTLPPGGAEAVDEYLEYRRIVGEDDGGKLFTPEEYEEYKKKVLPMRLQNRLFVSWRSPTGMDCKLVG PETLCFCTHRYKQHKTDFETIPQQRPIALPCRVSGCRCKAYHYVPLNGTQPIRCRCKHFADQHSAAALGFT CNACSKCSGFHSCYTCACGQPAYAHDTVETRQERLAQGKPVGRDVPYAAMGGTLGFSSLAEGYMRLLD S GIGAPSVEVLDSAVSAMDHPFLRAMHAPSTSSPQPLAGGNEDKAGKGC  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	28.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001165687</a></u>
Locus ID:	231946
UniProt ID:	<u><a href="#">Q8C790</a></u>


[View online »](#)

RefSeq Size:	2013
Cytogenetics:	6 B2.3
RefSeq ORF:	774
Synonyms:	D330028D13Rik