

## Product datasheet for **TP503401**

### Hsd17b10 (NM\_016763) Mouse Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Purified recombinant protein of Mouse hydroxysteroid (17-beta) dehydrogenase 10 (Hsd17b10), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug

**Species:** Mouse

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >MR203401 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MAAAVRSVKGLVAVVTGGASGLGLATAKRLVGQGATAVLLDVPDSEGEAQAKKLGESCIFAPANVTSEKE  
IQAALTLAKEKFGRIDVAVNCAGIAVAIKTYHQKKNKIHTLEDFQRVINVNLIQTFNVIRLVAGEMGQNE  
PDQGGQRGVIINTASVAAFEGQVGQAAYSASKGGIVGMLPIARDLAPTGRVVTIAPGLFATPLLTLPL  
EKVRNFLASQVPFPSRLGDPAEYAHVQTIENPFLNGEVIRLDGAIRMQP

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-MYC/DDK

**Predicted MW:** 27.3 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C after receiving vials.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_058043](#)

**Locus ID:** 15108

**UniProt ID:** [Q99N15](#)

**RefSeq Size:** 938



[View online »](#)

Cytogenetics: X F3

RefSeq ORF: 786

Synonyms: 17bHSD10; Ads9; ERAB; Hadh2