

Product datasheet for **TP503440**

Orc6 (NM_001163791) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse origin recognition complex, subunit 6 (Orc6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203440 protein sequence Red =Cloning site Green =Tags(s)
	<p>MESELVRRRLAPRLGLAEPVLRKAEFLRLSKVKCVLSARSSETSNAVICDLAASCRKCPDRAYLIR LSGLNKMVYQSLKSFECLLGLNSNVGIRDLAVQFSCTEAVNLAAEILQSYESGLPETQRADLDSRPLF TTAALLSACKILKIKVDKTKMITASGVKKAILDRLCKQLEKIGQQINRDSADLARPALKRKKPEFSPTLK KKEPGLPEPPAKEIEVIETLHKLPKDEDLTQDYEEWKRKILENAAKAQTATAE</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	29.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001157263</u>
Locus ID:	56452
UniProt ID:	<u>Q3TQX1</u>
RefSeq Size:	1533



[View online »](#)

Cytogenetics: 8 C3

RefSeq ORF: 789

Synonyms: 6720420110Rik; Orc6l