

Product datasheet for **TP503486**

Tmem41a (NM_025693) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse transmembrane protein 41a (Tmem41a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203486 protein sequence Red =Cloning site Green =Tags(s)
	 MRALLGLLLVFGGCTFALYLLSTRPLGPRLAAAGEPEGRSLWFPSDLAELRELVREYRKEHQAYVFLFCSAYLYKQGFAIPGSSFLNVLGALFGPWLGLLLCCVLTSVGATGCYLLSSLFGKQLVISYFPDKVALLQKKVEENRNSLFFLLFLRFPMPNWFNLNSAPILNIPVQFFSVLIGLIPYNFICVQTGSILSTLTSLDALFSWETVLKLLAIALVALVPGTLIKKFSQKRLALSETSDIGHPDRRKDA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	29.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_079969</u>
Locus ID:	66664
UniProt ID:	<u>Q9D8U2</u> , <u>Q9CUY8</u>
RefSeq Size:	1359



[View online »](#)

Cytogenetics: 16 B1

RefSeq ORF: 795

Synonyms: 2900010K02Rik; 5730578N08Rik