

Product datasheet for **TP503491**

C1qtnf6 (NM_028331) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse C1q and tumor necrosis factor related protein 6 (C1qtnf6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203491 protein sequence Red =Cloning site Green =Tags(s)
	 MRVIMGIASLGFLWAVFLLPLVFGVPTEETTFGESVASHLPKGCRRCCDPEDLMSSDDTVQAPVSPYVLP EVRPYINITILKGDKGDRGPTGTPGKPGKNGTRGDRGSQGVKGDKGQAGSPGSSCQTHYSAFSVGRKTGL HSENFLSLLFDRVFNVDGHFDMATGSFVAPLRGLYFFSLNVHSWNYKETYVHIVHNEQAVVILYAQPS ERSIMQSQSVMLPLVPGDRVWVRLFKRERENGIYSDVDVTYITFSGHLIKAEDN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	29.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_082607
Locus ID:	72709
UniProt ID:	Q6IR41 , A0A3B0IP29
RefSeq Size:	2443



[View online »](#)

Cytogenetics: 15 E1

RefSeq ORF: 795

Synonyms: 2810036M19Rik; CTRP6