

Product datasheet for **TP503587**

Rpa2 (NM_011284) Mouse Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Mouse replication protein A2 (Rpa2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug |
| Species: | Mouse |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >MR203587 protein sequence Red =Cloning site Green =Tags(s) |

MWNSGFESFSSSTYGGAGGYTQSPGGFGSPTPSQAEKKSrvRAQHIVPCTISQLLSATLTDEVFRIGDVE
ISQVTIVGIIRHAEKAPTNIYKIDDMTAPPMDVRQWVDTDDASGENAVVPPETYVKVAGHLRSFQNKKS
LVAFKIIPLEDMNEFTAHILEVNVSHMMLSKSNSQASAGRPSMSNPGMSEPGNFSGNFMFANGLTVVQN
QVLNLIKACPRPEGLNFQDLRSQLQHMPVPSIKQAVDFLCNEGHIYSTVDDDDHFKSTDAE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|----------------|--|
| Tag: | C-MYC/DDK |
| Predicted MW: | 29.4 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C after receiving vials. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_035414 |
| Locus ID: | 19891 |
| UniProt ID: | Q3TE40 |
| RefSeq Size: | 1782 |



[View online »](#)

Cytogenetics: 4 D2.3

RefSeq ORF: 813

Synonyms: AA409079; AI325195; AU020965; Rf-A2; RPA34