

## Product datasheet for TP503640

### Hccs (NM\_008222) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse holo cytochrome c synthetase (Hccs), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203640 representing NM_008222 Red=Cloning site Green=Tags(s)

MGASASSPATAVNASNASDGGQPASPPSGCPMHKGQRKGCPTAATSDLTSESKAHTVPAHQDRAYDYVEC  
PVTGARAKDKESLDPSNLMPPPNQTPSPDQPFTLSTSREESSIPRADSEKKWVYPSEQMFWNAMLRKGWK  
WKDDDISQKDMYNIIRIHNQNNEQAWKEILKWEALHAHECPCGSLVRFGGKAREYSPRARIRSWMGYEL  
PFDRHDWIINRCGTEVRYVIDYYDGGEVNKEYQFTILDVRPAFDSFSAVWDRMKVAWWRWTS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	31 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_032248</a>
Locus ID:	15159
UniProt ID:	<a href="#">P53702</a> , <a href="#">Q8BP79</a>
RefSeq Size:	2383



[View online »](#)

Cytogenetics: X F5

RefSeq ORF: 816

Summary: Links covalently the heme group to the apoprotein of cytochrome c.[UniProtKB/Swiss-Prot Function]