

Product datasheet for TP503695

Odf3l1 (NM_198673) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse outer dense fiber of sperm tails 3-like 1 (Odf3l1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203695 protein sequence Red =Cloning site Green =Tags(s) MKQSKGSRNYFYAQHPKEKEEVPSWHEIKQTPVIMATIKGPGPAKYLRSSCTGYIAHDISMFQEPAYSL HTRHTKKRIIDINSPGPCYFLNPKVTRFGISTCPQVPMEEERISNPRINCMASCKYNLEKTRPSGERQPP QYTFGYRCPYRVMDPNPAPNQYQLPVTLTGNIPVFRAAPSYSLASTNKNWFHKENIAGGPGPAMHTRPE P SVYQNRSPFLSMAKRFGCPLDHTHRPGPGSHDIQPVTVHKPRIPAFMTGKIHSPHLCPLIVDICD TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	31.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_941075</u>
Locus ID:	382075
UniProt ID:	<u>Q810P2</u>


[View online »](#)

RefSeq Size:	1080
Cytogenetics:	9 B
RefSeq ORF:	825
Synonyms:	BC049697; Gm1116