

Product datasheet for **TP503700**

Necap1 (NM_026267) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse NECAP endocytosis associated 1 (Necap1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from Mouse cDNA ORF Clone, MR203700, encoding Mouse full-length Necap1.
Tag:	C-MYC/DDK
Predicted MW:	30.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_080543
Locus ID:	67602
UniProt ID:	Q9CR95 , Q0VB06
RefSeq Size:	2393
Cytogenetics:	6 F2
RefSeq ORF:	825
Synonyms:	1200016B17Rik; AI747569; mFLJ00061
Summary:	Involved in endocytosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)