

Product datasheet for TP503734

Ccdc117 (NM_134033) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse coiled-coil domain containing 117 (Ccdc117), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203734 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAALGRPFSGLPPLSGSADFLQPPPAFAGRAFPFGAAGHDLAPRPGVVRGAPSSPGGRTARGRVSIIHCRKKH KRLAEDDECPVRKKRLTEAELGAVTDEWALGAHQGREGHGVNTCPSSLSMPSMLDVVCEEMDQTTGEPQC EVARRRLQEIEDRIIDEDEEVESDRNVSHLPSLVLSDTMKTGLKREFDEVFTKRMIESMSRPSMELVLWK PLPELLPEKPKPSSSPKNYRRESQAKHAAPGTAFPQRTEGLLEPQCADAPLYRSLEAATSTEEEMEL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	30.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_598794
Locus ID:	104479
UniProt ID:	Q6PB51
RefSeq Size:	3418



[View online »](#)

Cytogenetics: 11 A1

RefSeq ORF: 834

Synonyms: 1110004K02Rik; 1700026O03Rik; AU018638; AU042822; AV173073; BC018601