

Product datasheet for TP503780

Vsig4 (NM_177789) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse V-set and immunoglobulin domain containing 4 (Vsig4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203780 protein sequence Red =Cloning site Green =Tags(s)
	<p>MEISSGLLFLGHLIVLTYGHPTLKTPESTGTWKGDVKIQCIYDPLRGYRQVLVVKWLVRHGSDSVTIFLR DSTGDHIQQAKYRGRLKVKVPGDVSLQINTLQMDDRNHYTCEVTWQTPDGNQVIRDKIILRVRKYNP PRINTEAPTTLHSSLEATTIMSSTDLTTNGTGKLEETIAGSGRNLPIFAIFIISLCCIVAVTIPYILF RCRTFQQEYVYGVS RVFARKTSNSEETTRVTTIATDEPDSQALISDYDDPCLSQEYQITIRSTMSIPAC</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	31.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_808457
Locus ID:	278180
UniProt ID:	F6TUL9
RefSeq Size:	1432



[View online »](#)

Cytogenetics: X C3

RefSeq ORF: 843

Synonyms: A530061A11; BC025105; CRlg; Z39IG