

Product datasheet for TP503977

Oser1 (NM_025699) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse oxidative stress responsive serine rich 1 (Oser1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR203977 protein sequence Red =Cloning site Green =Tags(s) MKSEAKDGEEESLQTAFKKLRVDASGSIISLVGEGPSVRASARTAADDTKPKTMCASKDSWHGSTRKSS RGAVRTQRRRRSKSPVLHPPKFIHCSTTAPPSSSQLKHSQTEPPDGISGRGISTPKEFNAGENSTSLDV NHTGAAIEPLRSSVLRLPSESKTEELSDATQVSQESLTANDLSDFQSVSKLSQGKPCVCVGKECQCKRWH DMEVYSFSGQLQNPPLAPERRSLEDYSQSLHTRTLGSPRSCSEQARVYVDDVTIEDLAGYMEYYLYIPK KMSHMAEMMYT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	32 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_079975</u>
Locus ID:	66680
UniProt ID:	<u>Q9D722</u>


[View online »](#)

RefSeq Size:	1583
Cytogenetics:	2 H3
RefSeq ORF:	873
Synonyms:	3230401D17Rik; AI256769; Perit1