

## Product datasheet for **TP504138**

### Ogn (NM\_008760) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse osteoglycin (Ogn), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR204138 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 METVHSTFLLLLFVPLTQQAPQSQLDSHVNYEYATGNSEETKFSQDYEDKYLDGKSIKEKETMIIPDEKS LQLQKDEVIPSLPTKKENDEMPCLLCVCLSGSVYCEEVDIDAVPPLPKESAYLYARFNKIKKLTAKDFA DMPNLRRLDFTGNLIEDIEDGTFSKLSLLEELTLAENQLLRPLVLPKLTLLNAKHNKIKSGIKANTFK KLNKLSFLYLDHNDLESVPPNLPESLRVIHLQFNSSSLTDDTFCKANDTRYIRERIEEIRLEGNPIALG KHPNSFICKRLPIGSYF  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	34 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_032786</a>
Locus ID:	18295
UniProt ID:	<a href="#">Q62000</a> , <a href="#">Q543C5</a>



[View online »](#)

<b>RefSeq Size:</b>	2540
<b>Cytogenetics:</b>	13 A5
<b>RefSeq ORF:</b>	897
<b>Synonyms:</b>	3110079A16Rik; mimecan; mimican; OG; OIF; SLRR3A
<b>Summary:</b>	Induces bone formation in conjunction with TGF-beta-1 or TGF-beta-2.[UniProtKB/Swiss-Prot Function]