

Product datasheet for **TP504310**

Armc10 (NM_026034) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse armadillo repeat containing 10 (Armc10), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR204310 protein sequence Red =Cloning site Green =Tags(s)
	<p>MGGARDVGVWAAGLVLGAGACYCIYRLTRGPRRGRRLRPSRSAEDLTDGSYDDILNAEQLKLLYLLES TDDPVITEKALVTLGNNAFSTNQAIIRELGGIPIVGNKINSLNQSIEKALNALNNLSVNVENQTKIKI YVPQVCEDVFADPLNSAVQLAGLRLLTNMTVTNDYQHLLSGSVAGLFHLLLLGNGSTKVQVLKLLNLSE NPAMTEGLLSVQVDSSFLSLYDQGQVANEILLRALTLFQNIINCLKVEGRLANQIPFAKGSLLFFLYGEEC AQKMRALACHHDVDVKEKALAIKPKF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	33.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_080310
Locus ID:	67211
UniProt ID:	Q9D0L7



[View online »](#)

RefSeq Size:	2014
Cytogenetics:	5 A3
RefSeq ORF:	921
Synonyms:	2810037C14Rik; Svh
Summary:	May play a role in cell survival and cell growth. May suppress the transcriptional activity of p53/TP53 (By similarity).[UniProtKB/Swiss-Prot Function]