

Product datasheet for TP504437

Yipf2 (NM_138303) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse Yip1 domain family, member 2 (Yipf2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR204437 protein sequence Red =Cloning site Green =Tags(s)

MAAADELAFHEFEEATNLLAETPDAATTSQSDELTSREHVAVVWGSGIGYGAEVGEEEDDKTSLQEEKP
QPRFWTFDYQSFVDVTSQVLDRIKGSLLPHPGHNFVRHHLRNRPDLYGPFWICATLAFVLAVTGNLTL
VLAQRDPDSIHYSQFHKVTIAGITIYCYAWLVLALWGLRWRQGTRERMGLYTFLETVCVYGYSLFV
IPTVWLWLPVQWIQWLFGLGLALSAAGLVFTLWPVVREDTRLVAAALLSTVLLHALLALGCKLYFFQ
PLPLDHVVPAPQATPPSPNVLLPSSIQPMTTT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	34.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_612176</u>
Locus ID:	74766
UniProt ID:	<u>Q99LP8</u>



[View online »](#)

RefSeq Size: 2177
Cytogenetics: 9 A3
RefSeq ORF: 939
Synonyms: 1300010K09Rik